


PERFORMANCE

HS/HSG Catheter

OB-GYN



- Used for Hysterosalpingography and Sonohysterography procedures
- THIS PRODUCT DOES NOT CONTAIN LATEX
- Sterile. Disposable. Ready to use.
- Placement sheath decreases the need for a tenaculum
- Patient ambulatory once instrument properly inserted
- Allows for clear view of external cervical os
- Stylet available for cases when additional stiffness is needed

 **PERFORMANCE**Medical, Inc.

2512 Northland Drive | St. Paul, Minnesota 55120

www.Performance-MasterMedical.com

651.698.1631 | 800.798.1631

651.698.2083 fax